


OPTI QUIP

1200 SERIES POWER SUPPLIES

- RUNS 3 POPULAR LAMPS:
--75 WATT XENON
--100 WATT MERCURY
--200 WATT MERCURY
- LOW RIPPLE (5%)
- AUTOMATIC IGNITION
- COMBINED WATT-AMP-VOLT METER
- SOCKET FLEXIBILITY
Accommodates many popular lamphouses
- EXTREMELY QUIET
- LAMPHOUSE SAFETY INTERLOCK (MODEL 1205)
- RE-IGNITE HOT MERCURY BULBS IN 20 SECONDS



The features built into the Opti Quip 1200 DC power supplies are the ones that discriminating microscope users have demanded for over 25 years.


The 1200 allows the user to set and maintain the correct power for the lamp installed. Aging and the eventual need for replacement can be determined by observing the voltage at which the lamp is operating. A new lamp will operate at about 20 volts and will need replacing at about 35 volts.

The lamphouse safety interlock (Model 1205) allows operation only when the

socket is correctly installed and properly grounded to the power supply. Any attempt to open the lamphouse or to remove the socket while the lamp is operating will shut down the supply.

Users desiring newer technology and the ability to use a wider range of lamps (nine in all) should inquire about the Opti Quip 1600 series of power supplies. The 1600 series will also drive the Opti Quip Long Term Stabilizer (Model 1962).

See the back side of this sheet for technical specs and socket availability.



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

OPTI QUIP

1200 SERIES POWER SUPPLIES

Lamps powered:

75 watt xenon
 100 watt mercury
 200 watt mercury (but not the BA-027 lamp with 1 mm. arc)

Lamphouses driven:

The 1200 series will provide power for the Opti Quip Model 770 lamphouse, but sockets are available so that a variety of other lamphouses made by microscope manufacturers may be used. When ordering, specify the lamphouse that will be used so that an appropriate socket may be purchased also.

| Opti Quip Part Number | Lamphouse Manufacturer | Lamphouse Model | Lamp Type | Watts |
|-----------------------|------------------------|-----------------|-----------|-------|
| OQ315 | Leitz | 102 | xenon | 75 |
| universal | | 102 | mercury | 100 |
| OQ160 | | 100Z | mercury | 100 |
| OQ111* | | 250 | mercury | 200 |
| OQ526** | | 250 | xenon | 75 |
| OQ235Z | Zeiss | 100Z | xenon | 75 |
| OQ260 | | 100Z | mercury | 100 |
| OQ175* | | Universal | mercury | 200 |
| OQ530** | | Universal | xenon | 75 |

Adaptations:

*OQ125 Allows use of 100 watt mercury lamp in any socket that accepts a 200 watt mercury lamp
 **OQ513 Allows use of 75 watt xenon lamp in any socket that accepts a 150 watt xenon lamp
 OQ1680 Adaptation of Nikon HMX lamp socket (part #87531) to accept HBO 100W/2 or XBO 75W/2 (safety interlock version is optional)
 OQ1560 Modification of BX mercury arc lamp socket to accept either HBO 100W/2 mercury lamp or XBO 75W/2 for use with Opti Quip power supplies

Specifications:

Power required: 117V \pm 10% 50/60 Hz
 Power output: up to 200 watts
 Open circuit voltage: up to 75V D.C., dependent on lamp
 Maximum continuous current: 5.4 amps w/ 75 watt xenon Fuse: 5 amp
 Current ripple: 5% R.M.S.
 Line regulation: \pm 10% input causes $< \frac{1}{2}$ % output
 Ignition: automatic (Manual ignition is optional; see ATO on price list.)
 Trigger: 22K volts
 Size: 6" x 11" x 9"
 Weight: 11 pounds



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

Note: Table of Contents next page.

Thank you for purchasing an Opti-Quip 1500-1200 series power supply.

This manual is for use with the following supplies:

MODEL 1500 POWER SUPPLY

This supply will run any one of six different D.C. arc lamps - both xenon and mercury.

MODEL 1520 POWER SUPPLY

This supply is a variation on the standard 1500 which includes an extra receptacle with connections to power a 110V A.C. lamphouse mounted cooling fan and connections for lamphouse safety interlock. Any attempt to open the lamphouse makes the power supply automatically inoperative.

MODEL 1505 POWER SUPPLY

This supply has connections in its standard receptacle for a lamphouse safety interlock. Any attempt to open a lamphouse or remove a socket from a lamphouse designed with a safety interlock makes the power supply automatically inoperative.

MODEL 1200 POWER SUPPLY

This supply will run three different D.C. arc lamps - 75 watt xenon, 100 watt mercury and 200 watt mercury.

MODEL 1205 POWER SUPPLY

Same modifications as 1505 power supply.



TABLE OF CONTENTS

| | | |
|------------|--------------------------------------|--------------------|
| I | <u>FEATURES</u> | <u>PAGE</u> |
| | List of Features | 3 |
| | Front Panel Organization | 4 |
| | Terminology | 5 |
| | Volt-Amp-Watt Meter | 6 |
| | Hour Meter | 6 |
| | Fuse | 7 |
| II | <u>INSTALLATION</u> | |
| | Initial set-up | 8 |
| | Socket and Bulb | 9-10 |
| | Reading Volt-Amp-Watt Meter | 10 |
| | Bulb Selector Switch | 10 |
| | Hi-Low Trigger Switch | 11 |
| | "How to Adjust" | 12 |
| III | <u>OPERATION</u> | |
| | 200 Watt Mercury | 13 |
| | 200 Watt Mercury Stabilization | 14 |
| | Bulb Aging | 15 |
| IV | <u>PROBLEMS AND SERVICING</u> | |
| (A) | If servicing is necessary | 16 |
| (B) | Possible Problems | 16 |
| | Arcing over | 16 |
| | Blowing Fuses | 17 |
| | Intermittent Operation | 17 |
| | R.F. | 18 |
| | Problems & Service at a Glance | 19-20 |
| | A Final Note | 21 |
| (C) | Technical Assistance | 21 |
| | Integrated Circuit Description | 21 |
| | Schematics | 23-26 |
| | Parts List | 27-31 |



FEATURES

Introduction:

The Opti-Quip Solid State D.C. Power supplies have been designed to provide the finest illumination available for fluorescence and general microscopy.


Features of the Model 1500-1200 series D.C. Power Supplies

1. The model 1500 will run any one of six different arc lamps - 50-100-200 watt D.C. mercury lamps and 75-150-150W/GS (short gap) watt xenon lamps.
2. The 1200 can be used with a 100 watt mercury bulb, 75 watt xenon bulb or a 200 watt mercury bulb.
3. The life of the 200 watt D.C. mercury lamp is double that of the A.C. lamp.
4. The D.C. eliminates bulb flicker and gives more light and life than similar A.C. lamps.
5. The watt meter allows the user to adjust the power supply correctly for the particular bulb in use, and therefore maximize bulb life expectancy.
6. Fast Warm-up. Mercury lamps come up to full power in four minutes; xenon lamps in two minutes.
7. A non-resettable meter logs hours of use. This is helpful in establishing warranted bulb usage.
8. A special regulating circuit insures that $\pm 10\%$ variations in line voltage result in less than a 1% change in output power.
9. Power supplies are designed to re-ignite hot mercury bulbs in twenty seconds.
10. Opti-Quip manufactured sockets are compatible with a large number of existing microscope lamphouses.
11. The volt meter allows the user to monitor the aging process of the bulb and therefore plan for timely replacement.



FEATURES

1500 - 1200 Front Panel Organization



1. Pilot Light (lights up when supply is turned on)
2. Power switch (up-on position: down-off)
3. Fuse holder
4. Power adjust
5. Hour meter
6. Watt meter
7. Upper scale (reads watts and indicates correct wattage for mercury lamps)
8. Lower scale (reads amps or volts and indicates correct amps for xenon lamps)
9. Switch position to read "amps"
10. Switch to read "volts"
11. Switch to read "watts"
12. Switch
13. Bulb selector switch (1500 only)



FEATURES

Terminology:

- | | |
|------------------|--|
| Current | - Amps |
| Voltage | - Potential |
| I.C. | - Integrated Circuit |
| Sola Transformer | - Device used to maintain constant line voltage. |
| Plasma Streaming | - A magnetic field formed when mercury gas atoms ionize causing a cloud of positively charged atoms to form in a characteristic umbrella shape. |
| Arcing Over | - Trigger leaking. This may occur when the power supply is triggering. |
| Arcing | - The misdirected flow of current. |
| Power Adjust | - Recessed slotted device which is located at the front, lower left panel of the power supply. It can be rotated to the left or right with a screwdriver. In this way the watt or amp readings can be lowered, raised or set precisely on the meter. |
| Watts | - Power supplied to bulb - arrived at by multiplying volts times amps. |



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060


E-mail : info@meijitechno.com

<http://www.meijitechno.com>

FEATURES

Volt-Amp-Watt Meter: Mercury lamps are rated to be run at constant wattage. The operating voltage increases during the life of the bulb and necessitates a concurrent reduction in operating current. The watt meter integrates these two variables and allows the correct wattage to be set by use of the power adjust (see #4 page 4). As stated in the "High Pressure Xenon Lamp" catalog, "The maximum current for each xenon lamp must never be exceeded as otherwise safety of operation is no longer warranted, and the lamp life will be reduced essentially". If for some reason the watt meter reading doesn't register on the scale of your power supply, you can use the volt and amp meter readings to see what the watt reading would be. (Volts times amps equals watts). Osram also states that "5% either way (lamp running at too high or too low wattage) decreases the life of the bulb". With practice it also is possible to gauge the age of the particular bulb in use, with the aid of the Volt meter - Amp meter - Watt meter. A new 100 watt mercury lamp runs at about 20 volts. As the bulb ages, the arc length grows and the power supply must supply a higher voltage to maintain the arc. While maintaining 100 watts the operating current therefore drops. Older bulbs may go to as much as 26 amps at 38 volts. It is this lengthening of arc that is the real aging process.

Hour Meter: All Opti-Quip Power Supplies are equipped with a non resettable hour meter. This meter is connected for a portion of the testing before shipment. Right digit reads in tenths of an hour. The hour meter is a help to the customer to verify the warranted hours of bulb life. It is of importance to Opti-Quip if repairs become necessary, to have an idea of the total usage the power supply has received.



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA


Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

FEATURES



Fuse: The power supply main line is fused with a front panel fuse holder. In the event it is necessary to change the fuse, push in and down on the fuse holder to remove it from the supply. Fuse can then be replaced in the holder. Replace only with MTH type fuses (formerly 3AG) 8 amp for model 1500. 5 amp for model 1200.



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

INSTALLATION

PLEASE READ THE ENTIRE MANUAL CAREFULLY BEFORE USING THE POWER SUPPLY

This will insure perfect operation of your supply.

Supply Installation: Initial Set-up:

1. Plug the power supply into a 110V grounded outlet.
2. Set the dial below the watt meter to Volts position.
3. Turn on the power switch. Pilot light should come on. If volts register, you know the supply is getting power and will run the burner once it has triggered. (The meter should read 60-75 volts prior to ignition). (If pilot light does not light, check the fuse and replace it if necessary. See "Problems" page 16)

WARNING

Model 1200 uses 5A fuse, Model 1500 uses 8A fuse. Don't use a larger size fuse under any circumstance.

Triggering: Approximately 5 to 15 seconds after turning on the power switch, the just audible sound of "ticking" will be heard. This is a High Voltage trigger which ignites the bulb. TURN OFF THE POWER SWITCH. You are now assured that your power supply is capable of supplying power and should work for you. UNPLUG THE SUPPLY FROM THE OUTLET before installing the bulb in the socket and proceed.



INSTALLATION

Socket and Bulb Installation:

1. Install the bulb correctly in your socket. Reversed polarity can dramatically shorten the life expectancy of the particular bulb to as little as one or two minutes. The writing on the base of the bulb goes down (fits into the bottom of the socket) on both the mercury and xenon lamps.

Please note the Model #315 Model #406 Model #614 Model #645 Model #680 lamp sockets are universal sockets. They will accept either the xenon 75 watt bulb or the mercury 100 watt bulb. This is possible because either bulb may be run horizontally. When installing either bulb in these sockets the fat, positive (+) end of the bulb goes into the clamp and is tightened with the allen wrench.

2. If you have any difficulty inserting the lamp into the socket, it may be that the stamping (lettering) on a bulb is too heavy. To insure getting good contact and making sure that the bulb fits down into the socket, a small file can be used to file down the stamping on the lamp. The lamp is inserted correctly if the stamping (lettering) end is down.
3. Make sure that all connections are clean and tight.
4. Carefully insert the socket into the appropriate lamphouse. Be careful to follow special instructions included with particular socket you wish to use.
5. If the switch under the Watt meter is set to "amps" during ignition and warm up, no reading will be observed until the bulb actually is lit. It will be noted that initially mercury bulbs draw more than their nominally rated current.

| | nominal | nominal |
|------------|--------------|-------------|
| HBO 50W/3 | 2.2 amps | 22 volts |
| HBO 100W/2 | 5.0 amps | 20 volts |
| HBO 200W/2 | 3.1-4.2 amps | 47-65 volts |

This is normal operation. Once bulb has run for 4 or 5 minutes use "How to Adjust" (page 12) to set correct operating wattage. Xenon lamps however, come up to their operation current almost immediately and should not exceed their ratings.



INSTALLATION

Socket and Bulb installation cont'd:

6. Plug the socket into the power supply.
7. Plug the power supply into the outlet.
8. Review "How to Adjust" instructions (page 12).

To Read Meter

- a) To Read Watts The top line of the scale is used to read watts. Set switch to "Watts" position. The top scale reads from 0 to 250. See "How to Adjust" (page 12) for correct reading.
- b) To Read Volts The bottom line of the scale is used to read either Volts 0-100 or Amps 0-10. When switch is set to "Volts" position, add 0 to indicated number for correct voltage. For example, a "2" really is 20 volts.
- c) To Read Amps Each major division on scale is 1 amp, each small division is 0.2. For example, the line for "XBO 75" is 5.4 amps (see page 12 "How to Adjust")

Bulb Selector Switch (1500 series only)

See page 4 - "Front Panel Organization" #13
Before operation, set bulb selector switch so setting indicates bulb to be used. The 150 watt position is used both for the standard 150 watt xenon and the short gap 150 watt xenon, only the current ratings differ for these bulbs. Be sure you know which 150 watt bulb you are using. If, inadvertently, the switch is set for the wrong lamp after start up, simply switch to correct position as soon as noticed. Warning - some missettings will result in no bulb ignition, others could result in drastic bulb life reduction or possible violent failure if left uncorrected.



INSTALLATION

Hi-Low Trigger Switch

The 1500 series of power supplies has a "Hi" - "Low" switch setting for trigger voltage.


When shipped this switch is in the high position. This setting is appropriate where xenon bulbs primarily are to be used.

Under certain circumstances where a 50 watt mercury bulb is to be used, this "Hi" setting may result in hard starting. In this case the "Low" setting should be used.

The switch position is shown on the parts layout page. It is on the extreme left side of the board. This may be reached after removing the cover from the 1500 series power supply. Please contact us before changing this switch.



OPTI-QUIP 1500 - 1200 POWER SUPPLIES



MEIJI TECHNO AMERICA
 5895 Rue Ferrari, San Jose, CA 95138, USA
 Tel : 408-226-3454 Fax : 408-226-0900
 Toll Free : 1-800-832-0060
 E-mail : info@meijitechno.com
 http://www.meijitechno.com

VOLTS 0-100
 AMPS 0-10
 WATTS 0-250

Allow mercury bulbs five minutes and xenon bulbs three minutes to warm up. Then use "power adjust" to obtain meter readings, as indicated on this chart.

NOTE:

* BULBS SO MARKED FOR USE WITH 1500 SERIES ONLY

OPERATION

1. If you haven't already done so read "How to Adjust" instructions. (see page 12)
2. Turn on the Power Switch.
3. The bulb may ignite immediately or may take several seconds to start. This depends on the age of the bulb and whether or not the bulb is hot or cool.

WARNING: IF AT THIS TIME, LOUD "SNAPPING" SOUNDS ARE HEARD FROM THE LAMPHOUSE, TURN OFF THE POWER SUPPLY IMMEDIATELY. THIS INDICATES A PROBLEM WHICH WILL BE DEALT WITH UNDER PROBLEMS AND SERVICING.

If the bulb seems to be going on and off repeatedly before igniting, you may take a screwdriver and raise or lower the power adjust. This will help the bulb to ignite more rapidly. It may be necessary to do this, especially with an older lamp. As the bulb warms up, the current will drop to normal operating range. This should take about 2 or 3 minutes. The bulb has ignited when amps or watts register on the watt meter.

4. Operation of D.C. Power Supply with 200 watt mercury lamp.

As the bulb approaches the normal operating temperature it will "stabilize". At this time the arc which normally starts from the tip of one electrode and the side of the other, will jump to a position between the tips of both electrodes. In order for this to happen, the power supply automatically raises the operating voltage while the current is dropping during warm-up. This voltage will continue to rise until the arc stabilizes. If operating correctly, as the voltage rises, the current will drop until the correct wattage is maintained.

PLEASE NOTE: Periodically check operating parameters while lamp is running "watts" for mercury lamps. Current "amps" for xenon lamps. If the reading is slightly off, you can use a screwdriver to turn the power adjust until you obtain the exact reading desired. See "How to Adjust" (page 12)



OPERATION

Bulb Stabilization:

* * **NOTE:** Plasma Streaming can be confused with flickering (** see page 5 - "Terminology").

An older 200 watt mercury bulb can show its age by failing to stabilize.

A. There are three methods that you can use to determine whether or not your bulb has stabilized.

1. Observe arc through appropriate filters. You can see if the arc is between the two electrodes. Never look directly at an operating arc bulb.
2. Watch the watt meter while bulb is warming up toward the 200 watt reading. The needle will drop down very slightly and then go up when the arc has been maintained.
3. Watch the reflection of light, if possible, on a wall or on ceiling. The light will go off and on very quickly when the bulb stabilizes.

The bulb may flicker repeatedly after it stabilizes. In this case, use the power adjust to adjust the watt reading up or down slightly and the flickering should stop. Difficulty with stabilization occurs occasionally with new bulbs, but usually disappears after several hours use. It is exceedingly rare with the 100 watt mercury bulb.

B. If the bulb doesn't stabilize:

1. If it is an old bulb, you may have some difficulty with stabilization. In this case, use the power adjust to help the bulb to stabilize.
2. If it is a new bulb, the bulb should run for two hours or so and will stabilize after that. The power adjust may be used in this case also to aid the bulb in stabilizing.
3. If bulb goes out, let it cool for one minute before re-igniting.
4. Bulbs must be run at the correct wattage.



OPERATION

Bulb Stabilization cont'd from page 14

If the bulb does not stabilize, the power supply will automatically raise the operating voltage and lower the operating current to a point where the arc will no longer be maintained. At this point the bulb will stop working and the power supply will automatically start triggering again. Note the amp (current) reading when the bulb does go off. Restart the bulb. Use the power adjust to assist an old bulb (with high voltage) to get started. The old bulb may go on and off while the supply is triggering. If you replace the bulb and there is no longer an on/off phenomenon, you will know that the bulb was indeed too old.

Bulb Aging: Initial start up ages a bulb faster than continual operation at rated power. Diminuation of light output and lack of stability are indicative of a need for replacement.

The operating voltage of a 100 watt mercury bulb is indicative of the actual aging process. New bulbs start life operating around 20 volts. Once the operating voltage has reached about 38 volts, they should be replaced.



PROBLEMS AND SERVICING

NOTE: The serial number of the power supply is on the back panel of the supply.
If the occasion should ever arise that your power supply needs servicing and must be returned to Opti-Quip, please note the serial number of your supply and return the socket that you use along with the power supply.

* * IMPORTANT: Power supply must be double boxed and packed with appropriate materials, to prevent shipping damage.

A. If Servicing Becomes Necessary:

1. Turn off the power supply.
2. Leave the supply plugged into the outlet for about five minutes after turning it off. This gives the capacitors a chance to discharge to neutral.
3. Unplug the supply from the outlet.
4. Disconnect the socket from the power supply.

B. Possible Problems:

1. "Arcing Over" (Trigger leaking)

When the supply has been turned on, if a loud noise is heard coming from the lamphouse, the supply should be turned off immediately. This "arcing over" may be caused by:

- (a) A defective socket (in this case the "snapping" sound will be heard coming from the lamphouse).
- (b) Haphazard socket installation may result in wires being too close to or touching the lamphouse. (Please refer to specific instruction sheets for particular socket) "Arcing over" will occur only when the supply is triggering. If this phenomenon occurs, turn the supply off immediately. If this happens excessively, an I.C. may burn out.



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

PROBLEMS AND SERVICING

B. Possible Problems cont'd:

2. Blowing fuses

The length of time it takes for a fuse to fail can be indicative of the reason for failure.

- (a) As soon as power switch is turned on - disconnect lamp socket, replace fuse. If problem goes away, look for short circuit between lamphouse and socket. If problem remains, socket and power supply should be returned for repair. If the power supply is new, examine packing material for indications of damage.
- (b) If failure occurs after power supply has started to trigger, proceed as in (a). Probable cause is in socket or its installation. Correct as appropriate.
- (c) If failure occurs after power supply has ignited bulb and bulb has warmed up, the probable cause is a loose connection. This will cause grossly erroneous readings in the watt meter. A very loose connection will cause arcing with rapid severe heating of adjacent parts. The top connection to bulb is the first place to look. Arcing here will cause metal parts to be blackened or turn dark blue. This over heating will greatly reduce the life of the bulb.

3. Intermittent operation

In a situation where heavy equipment such as a centrifuge or a cooler is connected to the same branch circuit as the power supply, it is possible that the lamp may momentarily go out when this type of equipment is being started. Where observed, the phenomenon may be eliminated by connecting the power supply to a Sola constant voltage transformer - one with a 400-500 watt rating should be used.



PROBLEMS AND SERVICING

B. Possible Problems cont'd:

4. R.F. - Radio Frequency Interference

Arc lamps emit R.F. during operation. This interference may affect closed circuit television pictures. This is characterized by horizontal interference lines on the video screen. Several things may be done to keep this to a minimum.

- (a) Where ever possible use a camera with its electronics in the same case.
- (b) Be sure interconnecting cable between camera and monitor is shielded (aluminum foil wrapped around cable is a help).
- (c) Re-position power supply and monitor to effect minimum interference.
- (d) Attempt to connect power supply to a circuit which is on the other side of the main from the closed circuit system -an electrician in your building will be able to figure this out for you. Development work is progressing in this area.

5. R.F. - During Triggering

Large amounts of R.F. are given off by the high voltage trigger used to ignite arc lamps. This R.F. goes away after the bulb is ignited. To minimize the possible problems in this area, we suggest:

- (a) During the initial start up, be sure lamp power supply is first item turned on, only then should auxiliary microscope equipment be turned on.
- (b) Follow all suggestions as in number 4, particularly if a computer is involved.





MEIJI TECHNO AMERICA
5895 Rue Ferrari, San Jose, CA 95138, USA
Tel : 408-226-3454 Fax : 408-226-0900
Toll Free : 1-800-832-0060
E-mail : info@meijitechno.com
http://www.meijitechno.com

SERVICE AT A GLANCE

| <u>SYMPTOM</u> | <u>SIGN</u> | <u>CORRECTION AND/OR CAUSE</u> | <u>SEE PAGE</u> |
|---|---|---|---------------------|
| Loud Noise is heard coming from the lamphouse | "arcing over" (trigger leaking) | defective socket wires too close to or touching lamphouse | 16 |
| Supply stops operating for no apparent reason | Blown fuse (no triggering sound heard from supply) Bulb hasn't stabilized (supply will be triggering) "Brown-out" (triggering will be heard from supply) Arcing (blows fuse) (no triggering); running high; can't bring power down | Replace with 5 amp or 8 amp use as appropriate Turn off for 5 minutes and restart. Use power adjust to encourage lamp to stabilize. Sola transformer Tighten loose connection | 7 12 17 17 |
| Blows fuse | fuse blows right away After 45 minutes or so | Short circuit in lamp housing. Try disconnecting lamp connector on rear panel and try again. If problem persists, it's probably a defective power SCR or diode. Loose connection on lamp or short circuit in housing due to thermal expansion. | 17 |

Lamp flashes but will not ignite.

Power adjust set too low.



MEIJI TECHNO AMERICA
5895 Rue Ferrari, San Jose, CA 95138, USA
Tel : 408-226-3454 Fax : 408-226-0900
Toll Free : 1-800-832-0060
E-mail : info@meijitechno.com
http://www.meijitechno.com

MEIJI TECHNO

SERVICE AT A GLANCE

| <u>SYMPTOM</u> | <u>SIGN</u> | <u>CORRECTION AND/OR CAUSE</u> | <u>SEE PAGE</u> |
|----------------|-------------|--------------------------------|-----------------|
|----------------|-------------|--------------------------------|-----------------|

| | | | |
|---|--|---------------------|--|
| Pilot light on No open circuit voltage No trigger voltage | | Replace IC3 RCA 334 | |
|---|--|---------------------|--|

| | | | |
|--|--|---------------------------------|----|
| Bulb doesn't stabilize (Mercury 200 watt) | | Bulb may be too old New bulb | 14 |
|--|--|---------------------------------|----|

| | | | |
|--------|-----------------------------|--|----|
| Arcing | misdirected flow of current | Loose wire; discoloring on nickel plating on bulbs; signs of welding | 17 |
|--------|-----------------------------|--|----|

Periodically check watt meter
reading to be sure that too
much current is not being
drawn, and tighten loose wires
before damage occurs.

* For detailed description

PROBLEMS AND SERVICING

B. Possible Problems cont'd:

6. A Final Note

Many times a conversation with our chief electrical engineer and the help of the following pages, a knowledgeable service person can quickly rectify a seemingly complex problem. Our engineer is available on Tuesday, Wednesday, and Thursday and our telephone number is (914) 928-2254.

C. Technical Assistance

1. Integrated Circuit Description

There are six I.C.'s used in the 1500 and 1200 power supplies which are easy to replace and can sometimes solve minor problems with the supply.

I.C. (1) #LM556 supplies the trigger pulses to the regulating SCR's.

I.C. (2) #LM324 drives I.C. (1), and performs all regulating functions, except for the current sense amplifier, which is in I.C. (3).

I.C. (3) #LM324 section (2) is the current sense amp. Sections (1) and (4) control the Hi voltage trigger, and section (3) prevents an over voltage condition prior to lamp ignition.

I.C. (4) #LM324 comprise the multiplying watt meter circuit

I.C. (6) #LM556

I.C. (5) #4016BCP is an electronic switch that selects the appropriate signal for the front panel meter.



PROBLEMS AND SERVICING

C. Technical Assistance cont'd:

2. When Replacing I.C.'s


- (a) Before removing the I.C. from the board, please note the polarity on the I.C. Polarity must be correct. Replace the I.C. in the same direction as you find it on the board.
- (b) Look for the mark on the I.C. (There may be more than one). This also indicates direction for replacement.
- (c) Get the prongs in straight. Do not bend them.

3. The Schematics - see successive four pages.

4. A Physical Parts Layout - see page 27.


5. A Parts List - see pages 28-31.





1200 WATTMETER
19 APRIL 84 MJD PAGE 2 of 3





OPT1 QUIP MODEL 1200

MEIJI TECHNO AMERICA
 5895 Rue Ferrari, San Jose, CA 95138, USA
 Tel : 408-226-3454 Fax : 408-226-0900
 Toll Free : 1-800-832-0060
 E-mail : info@meijitechno.com
<http://www.meijitechno.com>

MEIJI TECHNO

1200/1500 Series Electronics Parts List 02/04/87

| REFERENCE | 1200 | 1500 | 1520 | 1205 |
|-----------|--------|--------|--------|------------|
| R 1 | 10K | 10K | 10K | 10K |
| R 2 | 10K | 10K | 10K | 1/4W 1% |
| R 3 | .01 | .01 | .01 | 1/4W 1% |
| R 4 | 1Meg | 1Meg | 1Meg | 3W 5% |
| R 5 | 10K | | | 1/4W 1% |
| R 6 | 10K | | | 10K |
| R 7 | 50K | | | 1/4W 5% |
| R 8 | 30K | 30K | 30K | 50K |
| R 9 | 10K | 10K | 10K | TrimPot |
| R10 | 470K | 470K | 470K | 1/4W 5% |
| R11 | 150K | JUMPER | JUMPER | 10K |
| R12 | 100K | 100K | 100K | 1/4W 5% |
| R13 | 82K | JUMPER | JUMPER | 100K |
| R14 | 10K | 10K | 10K | POWER ADJ. |
| R15 | 2K | 1K | 1K | 1/4W 5% |
| R16 | 2.7Meg | 2.7Meg | 2.7Meg | 1/4W 5% |
| R17 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R18 | 10Meg | 10Meg | 10Meg | 1/4W 5% |
| R19 | 9.1Meg | 9.1Meg | 9.1Meg | 1/4W 5% |
| R20 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R21 | 100K | 100K | 100K | 1/4W 5% |
| R22 | 30K | 30K | 30K | 1/4W 5% |
| R23 | 10K | 10K | 10K | 1/4W 5% |
| R24 | 1K | 1K | 1K | 1/4W 5% |
| R25 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R26 | 2K | 2K | 2K | 1/4W 5% |
| R27 | 33.2K | 33.2K | 33.2K | 1/4W 1% |
| R28 | 3.01K | 3.01K | 3.01K | 1/4W 1% |
| R29 | 510K | 510K | 510K | 1/4W 5% |
| R30 | 10Meg | 10Meg | 10Meg | 1/4W 5% |
| R31 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R32 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R33 | 15K | 15K | 15K | 1/4W 5% |
| R34 | 100K | 100K | 100K | 1/4W 5% |
| R35 | 1K | 1K | 1K | 1/4W 5% |
| R36 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R37 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R38 | 180K | 100K | 100K | 1/4W 5% |
| R39 | 3.3K | 3.3K | 3.3K | 1/4W 5% |
| R40 | 5K | 5K | 5K | 1/4W 5% |
| R41 | 30K | 30K | 30K | TrimPot |
| R42 | 270K | 270K | 270K | 1/4W 5% |
| R43 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R44 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R45 | 330K | 330K | 330K | 1/4W 5% |
| R46 | 100K | 100K | 100K | 1/4W 5% |
| R47 | 1.5K | 1.5K | 1.5K | 1/4W 5% |
| R48 | 20K | 20K | 20K | 1/4W 5% |
| R49 | 50K | 50K | 50K | TrimPot |



| REFERENCE | 1200 | 1500 | 1520 | 1205 |
|-----------|----------|---------|---------|------------------|
| R50 | 2Meg | 2Meg | 2Meg | 1/4W 5% |
| R51 | 1Meg | 1Meg | 1Meg | 1/4W 5% |
| R52 | 604K | 604K | 604K | 1/4W 1% |
| R53 | 604K | 604K | 604K | 1/4W 1% |
| R54 | 604K | 604K | 604K | 1/4W 1% |
| R55 | 604K | 604K | 604K | 1/4W 1% |
| R56 | 47K | 47K | 47K | 1/4W 5% |
| R57 | 15K | 15K | 15K | 1/4W 5% |
| R58 | 5K | 5K | 5K | TrimPot |
| R59 | 5K | 5K | 5K | TrimPot |
| R60 | 6.8K | 6.8K | 6.8K | 1/4W 1% |
| R61 | 15K | 15K | 15K | 1/4W 5% |
| R62 | 15K | 15K | 15K | 1/4W 5% |
| R63 | 15K | 15K | 15K | 1/4W 5% |
| R64 | 10K | 10K | 10K | 10K 1/4W 5% |
| R65 | 510K | 510K | 510K | 510K 1/4W 5% |
| R66 | 100K | 100K | 100K | 100K 1/4W 5% |
| R67 | NOT USED | | | |
| R68 | 10K | 10K | 10K | 10K 1/4W 5% |
| R68 | 10K | 10K | 10K | 10K 1/4W 5% |
| R69 | 330 Ohm | 330 Ohm | 330 Ohm | 330 Ohm 1/2W 5% |
| R70 | 330 Ohm | 330 Ohm | 330 Ohm | 330 Ohm 1/2W 5% |
| R71 | 4.7K | 4.7K | 4.7K | 2 W 10% |
| R72 | 0.5 Ohm | 0.5 Ohm | 0.5 Ohm | 50W 10% |
| R73 | 1Meg | 1Meg | 1Meg | 1Meg 1/4W 5% |
| R74 | 1Meg | 1Meg | 1Meg | 1Meg 1/4W 5% |
| R75 | 4.3K | 4.3K | 4.3K | 4.3K 1/4W 5% |
| R76 | 33 Ohm | 33 Ohm | 33 Ohm | 33 Ohm 1/4W 5% |
| R77 | 4.3K | 4.3K | 4.3K | 4.3K 1/4W 5% |
| R78 | 15K | 15K | 15K | 15K 1/2W 10% |
| R79 | 15K | 15K | 15K | 15K 1/2W 10% |
| R80 | 100 Ohm | 100 Ohm | 100 Ohm | 100 Ohm 1/2W 10% |
| R81 | 10K | 10K | 10K | 10K 2 W 10% |
| R82 | 10K | 10K | 10K | 10K 2 W 10% |
| R83 | 3Meg | 3Meg | 3Meg | 3Meg 1/4W 5% |
| R84 | 6.8K | 6.8K | 6.8K | 1/4W 5% |
| R85 | 4.7K | 4.7K | 4.7K | 1/4W 5% |
| R107 | | 1K | 1K | 1/4W 5% |



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

| REFERENCE | 1200 | 1500 | 1520 | 1205 | |
|-----------|----------|----------|----------|----------|---------------|
| C 1 | 10ufd | 10ufd | 10ufd | 10ufd | 16V Tant. |
| C 2 | 10ufd | 10ufd | 10ufd | 10ufd | 16V Tant. |
| C 3 | 4.7ufd | 4.7ufd | 4.7ufd | 4.7ufd | 16V Tant. |
| C 4 | .0047ufd | .0047ufd | .0047ufd | .0047ufd | 100V Mylar |
| C 5 | 0.1ufd | 0.1ufd | 0.1ufd | 0.1ufd | 16V Tant. |
| C 6 | 0.1ufd | 0.1ufd | 0.1ufd | 0.1ufd | 16V Tant. |
| C 7 | .001ufd | .001ufd | .001ufd | .001ufd | 100V Mylar |
| C 8 | 4.7ufd | 4.7ufd | 4.7 | 4.7 | 16V Tant. |
| C 9 | .47ufd | .47ufd | .47ufd | | 100V Mylar |
| C10 | 0.1ufd | 0.1ufd | 0.1ufd | | 16V Tant. |
| C11 | 0.1ufd | 0.1ufd | 0.1ufd | | 16V Tant. |
| C12 | .022ufd | .022ufd | .022ufd | | 100V Mylar |
| C13 | .047ufd | .047ufd | .047ufd | | 100V Mylar |
| C14 | .047ufd | .047ufd | .047ufd | | 100V Mylar |
| C15 | 0.1ufd | | | | 16V Tant. |
| C16 | .01ufd | .01ufd | .01ufd | .01ufd | 1KV Ceram. |
| C17 | .01ufd | .01ufd | .01ufd | .01ufd | 1KV Ceram. |
| C18 | .01ufd | .01ufd | .01ufd | .01ufd | 1KV Ceram. |
| C19 | 17Kufd | 17Kufd | 17Kufd | | 75V Elect. |
| C20 | 0.1ufd | 0.1ufd | 0.1ufd | 0.1ufd | 600V Mylar |
| C21 | 5ufd | 5ufd | 5ufd | 5ufd | 150V Elect. |
| C22 | Jumper | 5ufd | 5ufd | Jumper | 150V Elect. |
| C23 | 15ufd | 15ufd | 15ufd | 15ufd | 16V Tant. |
| C24 | 4ufd | 4ufd | 4ufd | 4ufd | 350V Elect |
| C25 | 4ufd | 4ufd | 4ufd | 4ufd | 350V Elect. |
| C26 | .22ufd | .22ufd | .22ufd | .22ufd | 600V Polypro. |
| C27 | .22ufd | .22ufd | .22ufd | .22ufd | 600V Polypro. |
| C28 | 250ufd | 250ufd | 250ufd | 250ufd | 25V Elect. |
| C29 | 4.7ufd | 4.7ufd | 4.7ufd | | 16V Tant. |
| C30 | 4.7ufd | 4.7ufd | 4.7ufd | 4.7ufd | 16V Tant. |
| C31 | 0.1ufd | 0.1ufd | 0.1ufd | | 50V Ceram. |
| C32 | 0.1ufd | 0.1ufd | 0.1ufd | 0.1ufd | 100V Mylar |
| C33 | 0.1ufd | 0.1ufd | 0.1ufd | 0.1ufd | 50V Ceram. |
| C34 | 15ufd | 15ufd | 15ufd | 15ufd | 16V Tant. |
| C35 | .22ufd | .22ufd | .22ufd | .22ufd | 600V Polypro. |
| C36 | .22ufd | .22ufd | .22ufd | .22ufd | 600V Polypro. |
| C37 | 10ufd | 10ufd | 10ufd | 10ufd | 16V Tant. |
| C38 | .01ufd | .01ufd | .01ufd | .01ufd | 50V Mylar |



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

| REFERENCE | 1200 | 1500 | 1520 | 1205 |
|-----------|---------|-----------|-----------|-------------------------------|
| IC 1 | 556 | 556 | 556 | 556 |
| IC 2 | 324 | 324 | 324 | 324 |
| IC 3 | 324 | 324 | 324 | 324 |
| IC 4 | 324 | 324 | 324 | |
| IC 5 | 4066 | 4066 | 4066 | |
| IC 6 | 556 | 556 | 556 | |
| IC 7 | 7815 | 7815 | 7815 | 7815 |
| OC 1 | 4N26 | 4N26 | 4N26 | 4N26 |
| CR 1 | 1N914 | 1N914 | 1N914 | 1N914 |
| CR 2 | 1N914 | 1N914 | 1N914 | 1N914 |
| CR 3 | 1N4736A | 1N4736A | 1N4736A | 1N4736A 6.8Volt Zener |
| CR 4 | 1N4733A | 1N4733A | 1N4733A | 1N4733A 5Volt Zener |
| CR 5 | 1N914 | 1N914 | 1N914 | 1N914 |
| CR 6 | 1N914 | 1N914 | 1N914 | 1N914 |
| CR 7 | 1N914 | 1N914 | 1N914 | 1N914 |
| CR 8 | 1N4733A | 1N4733A | 1N4733A | |
| CR 9 | MR752 | 1/4 BR358 | 1/4 BR358 | Mounted on Chassis (BR358) |
| CR10 | MR752 | 1/4 BR358 | 1/4 BR358 | |
| CR11 | MR752 | MR752 | MR752 | MR752 |
| CR12 | 1N4007 | 1N4007 | 1N4007 | 1N4007 |
| CR13 | 1N4007 | 1N4007 | 1N4007 | 1N4007 |
| CR14 | 1N4007 | 1N4007 | 1N4007 | 1N4007 |
| CR15 | | 1N4007 | 1N4007 | |
| CR16 | 1N4007 | 1N4007 | 1N4007 | 1N4007 |
| CR17 | 1N4007 | 1N4007 | 1N4007 | 1N4007 |
| CR18 | 1N4745A | 1N4745A | 1N4745A | 1N4745A 15V Zener |
| CR19 | 1N4007 | 1N4007 | 1N4007 | 1N4007 |
| BR 1 | CSCB40S | CSCB40S | CSCB40S | CSCB40S 400V @ 1.5A |
| V 1 | V130LA1 | V130LA1 | V130LA1 | V130LA1 |
| V 2 | V130LA1 | V130LA1 | V130LA1 | V130LA1 |
| Q 1 | 2N5087 | 2N5087 | 2N5087 | 2N5087 |
| Q 2 | 2N3669 | 2N3669 | 2N3669 | 2N3669 RCA |
| Q 3 | 2N3669 | 2N3669 | 2N3669 | 2N3669 RCA |
| Q 4 | MPU131 | MPU131 | MPU131 | MPU131 |
| Q 5 | S2600M | S2600M | S2600M | S2600M RCA |
| Q 6 | 2N5087 | 2N5087 | 2N5087 | |
| T 1 | ST-3-16 | ST-3-16 | ST-3-16 | ST-3-16 |
| T 2 | TX107 | TX107 | TX107 | |
| L 1 | CH-4 | CH-6 | CH-6 | |
| L 2 | LN103 | LN103 | LN103 | LN103 |
| M 1 | 0-1Ma | 0-1Ma | 0-1Ma | |



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA


Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>

RELATIVE SPECTRAL DISTRIBUTION FOR MERCURY AND XENON



BIBLIOGRAPHY:

"Florescence Microscopy"; Ernst Leitz
 "Super Pressure Mercury"; Osram
 "XBO Xenon Short Arc Lamp"; Osram

— XENON CURVE

—*— MERCURY CURVE



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>


OPTI-QUIP 1500 1200 POWER SUPPLIES

| BULB | RATED TO RUN AT | METER SWITCH POSITION | METER READING |
|----------------------------|-----------------|-----------------------|---------------|
| 200 WATT MERCURY | 200 WATTS | - WATTS | |
| 100 WATT MERCURY | 100 WATTS | - WATTS | |
| * 50 WATT MERCURY | 50 WATTS | - WATTS | |
| 75 WATT XENON | 5.4 AMPS | - AMPS | |
| 100 WATT XENON | 6.7 AMPS | AMPS | |
| * 150 WATT XENON | 7.5 AMPS | - AMPS | |
| * 150 WATT XENON SHORT GAP | 8.5 AMPS | - AMPS | |

NOTE:


Allow mercury bulbs five minutes and xenon bulbs three minutes to warm up. Then use "power adjust" to obtain meter readings, as indicated on this chart.

* BULBS SO MARKED FOR USE WITH 1500 SERIES ONLY



MEIJI TECHNO AMERICA
 5895 Rue Ferrari, San Jose, CA 95138, USA
 Tel : 408-226-3454 Fax : 408-226-0900
 Toll Free : 1-800-832-0060
 E-mail : info@meijitechno.com
<http://www.meijitechno.com>

MEIJI TECHNO



MEIJI TECHNO AMERICA

5895 Rue Ferrari, San Jose, CA 95138, USA

Tel : 408-226-3454 Fax : 408-226-0900

Toll Free : 1-800-832-0060

E-mail : info@meijitechno.com

<http://www.meijitechno.com>