


Meiji Techno Boom Mounted Optics Part Numbers & Specifications



System# unimac1

Unlike stereo optics pods which are only using one of the optical paths to connect with a camera, the UNIMAC system is essentially one zoom optics element. And since UNIMAC is modular, there are hundreds of possible configurations, each one tailored to your imaging tasks.

Part Numbers & Specifications

Part Numbers		Description
A	S-4300	Universal Boom Stand requires FS block, vertical pole: 400mm, horizontal pole: 535mm. Base: 265mm x 265mm x 25mm
B	No Eyepieces	No eyepieces on this configuration, direct-to-camera
C	MS-3	Zoom objective, 0.7x to 4.5x
	MS-4	Intermediate Tube
	MS-5	Adapter Plate, accepts FS holder
	MS-1	Straight camera tube
D	FS	Inclinable and tiltable coarse focus block, accepts S-4300 and S-4400 boom stands.
E	FT190/115	150 watt Illuminator Power Supply with variable intensity control
	FL150/65	Mini-annular ring light guide
F	CC2100	USB2 Digital Video Camera, 1.4MP, complete with software



**MEIJI
TECHNO**

Meiji Techno America 3010 Olcott Street Santa Clara, CA 95054
Toll Free USA: 1(800) 832-0060 (408) 970-4799 FAX: (408) 970-5054
Visit us online at www.meijitechno.com