


Meiji Techno Boom Mounted Optics Part Numbers & Specifications



System# unimac3

This UNIMAC system incorporates an erect-image, trinocular viewing head with Cmount and digital video camera. The mini-annular ring light provides 150 watts of bright, even lighting for any imaging tasks.

Part Numbers & Specifications

Part Numbers		Description
A	CC2100	USB2 Digital Video Camera, 1.4MP, complete with software
B	MA151/10	Cmount adapter, no lens
C	MA407	10x eyepieces, fits MA655/05 Trinocular Erect Viewing Head
D	MA655/05	Trinocular Erect Viewing Head w/ beamsplitter & de-polarizer
E	MS-4	Intermediate Tube
	MS-5	Adapter Plate
	MS-3	Zoom objective, 0.7x to 4.5x
F	FS	Inclinable, tiltable coarse focus block accepts S-4300 and S-4400 boom stands
G	FL150/65	Mini-annular ring light guide
H	S-4300	Universal Boom Stand requires FS block, vertical pole: 400mm, horizontal pole: 535mm. Base: 265mm x 265mm x 25mm
I	FT190/115	150 watt Illuminator Power Supply with variable intensity control



Meiji Techno America 3010 Olcott Street Santa Clara, CA 95054
Toll Free USA: 1(800) 832-0060 (408) 970-4799 FAX: (408) 970-5054
Visit us online at www.meijitechno.com